Instruction Sheet

Human Airway Epithelial Cells

FC-0016	LL-0023
FC-0106	LM-0007
FC-0035	LS-1047
FC-0011	
FC-0054	
FC-0065	



This product is for Research Use Only. This product is not approved for human or veterinary use or for use in in vitro diagnostics or clinical procedures.





Better Solutions for Breakthrough Results FC-0016 LL-0023 FC-0106 LM-0007 FC-0035 LS-1047 FC-0011 FC-0054 FC-0065



Human Bronchial/Tracheal Epithelial Cells (HBTEC)	Human Small Airway Epithelial Cells (HSAEC)
Human Lobar Bronchial Epithelial Cells	Human Small Airway Epithelial Cells-Asthma +
(ILDEC)	(HSAEC-A+COPD)

Safety and Use Statement

Lifeline recommends storing cryopreserved vials in liquid nitrogen vapor phase. Handle cryopreserved vials with caution. Always wear eye protection and gloves when working with cell cultures. Aseptically vent any nitrogen from cryopreserved vials by carefully loosening the vial cap in a biosafety cabinet prior to thawing the vials in a water bath. If vials must be stored in liquid phase, the vials should be transferred to vapor phase storage or -80°C for up to 24 hours prior to being thawed.

Quick Steps for Cell Culture:

- 1. Always wash hands before and after working with cell cultures.
- 2. Always wear eye protection and gloves when working with cell cultures.
- 3. Store cryopreserved cells in liquid nitrogen vapor phase.
- 4. When working with cells or medium, always use a certified biological safety cabinet.
- 5. Handle cryopreserved vials with caution. Aseptically vent any nitrogen from cryopreserved vials in a biosafety cabinet prior to thawing in a water bath.
- 6. Feed cells using pre-warmed Lifeline® culture medium according to feeding chart.
- 7. When cells are 70 to 80% confluent and actively proliferating, passage cells using Lifeline subculture reagents.
- 8. Rinse cultures with PBS.
- 9. Trypsinize cells until rounded, do not over trypsinize.
- 10. Add TNS to stop trypsinization.
- 11. Add more PBS to rinse all the cells from the culture surface.
- 12. Centrifuge at 150 x g for 3 to 5 minutes. Adjust speed and time as appropriate for your centrifuge.
- 13. Aspirate solution from centrifuge tube, add pre-warmed Lifeline culture medium and gently re-suspend cells.
- 14. Count cells using a hemacytometer, re-plate at 5,000 cells per cm² (or use a split ratio of 1:4 to 1:10) in vessel containing pre-warmed Lifeline culture medium.



Call Lifeline Technical Service and Sales at 877.845.7787 Or visit lifelinecelltech.com for more information Lifeline Cell Technology | 8415 Progress Drive, Suite T | Frederick, MD 21701 ©2024 Lifeline Cell Technology, an International Stem Cell Company | All Rights Reserved



- 15. Incubate cells using 1 mL of culture medium per 5 cm² surface area at 37°C, 5% CO₂.
- 16. For any question on cell handling, please contact technical service. We are here to help.

Lifeline Technical Note: There are different and often contradictory terminologies used by cell culture companies to define the passage number of cells. Lifeline's designation of 'primary cells' are cells that have been isolated from tissue, plated onto culture vessels, expanded, harvested and cryopreserved. The term 'secondary' indicates that the cells have been isolated, plated and expanded in culture vessels twice before being harvested for cryopreservation.

Product		Catalog Number(s)
HSAEC	Normal Human Small Airway Epithelial Cells, Primary – 500,000 cells/vial	FC-0016
HSAEC	Normal Human Small Airway Epithelial Cells, Secondary – 500,000 cells/vial	<u>FC-0106</u>
HBTEC	Normal Human Bronchial/Tracheal Epithelial Cells, Primary – 500,000 cells/vial	<u>FC-0035</u>
HBTEC	Normal Human Bronchial/Tracheal Epithelial Cells, Secondary – 500,000 cells/vial	FC-0011
HLBEC	Normal Human Lobar Bronchial Epithelial Cells, Primary – 500,000 cells/vial	<u>FC-0054</u>
HSAEC- A+COPD	Diseased Human Small Airway Epithelial Cells-Asthma + Chronic Obstructive Pulmonary Disease. Primary – 500.000 cells/vial	FC-0067
BronchiaLife [*] (BronchiaLife	[™] Medium Complete Kit Basal Medium, BronchiaLife LifeFactors™ Kit)	<u>LL-0023</u>
BronchiaLife	™ Epithelial Basal Medium	<u>LM-0007</u>
BronchiaLife	[™] LifeFactors Kit	<u>LS-1047</u>
GA Antimicrobial Supplement, 0.5 mL (Gentamicin 30 mg/mL, Amphotericin B 15 µg/mL); provided with purchase of LL-0023		<u>LS-1104</u>
TrypKit™ Subculture Reagent Kit		<u>LL-0013</u>
Phenol Red Supplement [33 mM], 1 mL		<u>LS-1009</u>
FrostaLife™ C	ryopreservation Solution, 100 mL	LM-0015

Notes:

Basic Aseptic Technique

Cells, medium, and LifeFactors[™] should only be used in an aseptic environment, a Class II biological safety cabinet with front access and filtered laminar airflow, or an equivalent device. Always wear gloves and eye protection when working with these materials. Wipe or spray all bottles and vials with 70% ethanol or isopropanol, especially around the area of the cap, before placing them in the biological safety cabinet. Allow these surfaces to dry completely before opening the bottle or vials. Transfer cells, medium, or LifeFactors with disposable sterile pipettes. Do not mouth pipette! Take up the volume needed into the pipette, being careful not to touch the sterile tip to the rim of the container or any other surface. Close the container and open the container into which the transfer is being made, again being careful not to touch any surfaces with the sterile tip. Transfer the material and close the container. Wash your hands before and after working with cell cultures. Do not block airflow in a laminar flow hood as this may compromise sterility. Ensure that biological safety cabinets are certified routinely and that the HEPA filters are replaced regularly.



Call Lifeline Technical Service and Sales at 877.845.7787 Or visit lifelinecelltech.com for more information Lifeline Cell Technology | 8415 Progress Drive, Suite T | Frederick, MD 21701

Lifeline Cell Technology | 8415 Progress Drive, Suite T | Frederick, MD 21701 ©2024 Lifeline Cell Technology, an International Stem Cell Company | All Rights Reserved



Cryopreserved Vials

FC-0065

Human Airway Epithelial Cells are sold as cryopreserved vials and are shipped in biodegradable insulated packages containing dry ice to ensure the cells remain in a cryopreserved state. To maintain the cells' integrity, unpack the products immediately upon receipt and store at a temperature lower than -150°C or in the vapor phase of a liquid nitrogen dewar. If the cells are to be thawed and plated within 24 hours, they may be stored at -80°C. Do not store the vial for more than 24 hours at -80°C, as the cells will slowly degrade at this temperature.

Medium Preparation

BronchiaLife[™] Basal Medium (LM-0007) contains no growth factors, antimicrobials, or phenol red. To support Human Airway Epithelial Cell proliferation, you must add BronchiaLife LifeFactors (LS-1047) to the basal medium (see BronchiaLife Cell Culture Medium instructions).

A vial of Gentamicin and Amphotericin B (GA; LS-1104) is provided with the purchase of BronchiaLife Medium Complete Kit (LL-0023) for your convenience. The use of GA is recommended to inhibit potential fungal or bacterial contamination of eukaryotic cell cultures. Phenol Red (LS-1009) may be purchased, but is not required.

Pre-warming Medium

If using less than 100 mL of medium, Lifeline[®] recommends warming only the volume needed in a sterile conical tube. Repeated warming of the entire bottle over extended periods will cause degradation and reduced shelf life of the medium. Medium will warm to 37°C in 10 to 30 minutes, depending on the volume. Do not leave medium in water bath for extended periods.

Thawing and Plating Cryopreserved Cells

Remove vial from dewar and check the cap to be sure that the vial is securely sealed. <u>Spray the vial with 70%</u> ethanol or isopropanol and transfer it to a biosafety cabinet. Allow it to dry thoroughly and carefully loosen the cap to vent any liquid nitrogen that may have entered the vial. Recap the vial and hold only the bottom half of the vial in a 37°C water bath for approximately one minute or until vial is almost completely thawed—a small amount of ice should still be visible. To avoid potential contamination, do not allow the vial cap to make contact with the water. Do not over thaw as this may damage the cells. Dry the vial, spray the exterior of the vial with 70% ethanol or isopropanol and place the vial in a biological safety cabinet and allow it to dry. Carefully remove the cap to avoid contamination or splatter. Gently resuspend the cells in the vial using a 1- or 2-mL sterile pipette. Do not centrifuge; the cells may be directly plated from the vial. Plate cells into pre-warmed fully supplemented BronchiaLife[™] Medium in the desired culture vessel at a density of 5,000 cells per cm². Flasks with vented caps, commonly referred to as filter caps, are strongly recommended. Gently rock the culture vessel from side to side and front to back to evenly distribute cells within the vessel. Place seeded culture vessel in a 37°C, 5% CO₂ incubator. Refeed the cells after they have attached (approximately 4 to 36 hours after inoculation) to remove cryopreservation reagents.



Call Lifeline Technical Service and Sales at 877.845.7787 Or visit lifelinecelltech.com for more information Lifeline Cell Technology | 8415 Progress Drive, Suite T | Frederick, MD 21701 ©2024 Lifeline Cell Technology, an International Stem Cell Company | All Rights Reserved FC-0016 LL-0023 FC-0106 LM-0007 FC-0035 LS-1047 FC-0011 FC-0054

FC-0065



Standard Calculation for Plating of Cells

Gently re-suspend the cells evenly in the medium. Using a clean hemocytometer, and aseptic technique, remove 20 μ L of the cell suspension to a separate tube, such as a microcentrifuge tube. Add an equal volume of 0.4% Trypan Blue solution to the cell suspension in the microcentrifuge tube. Perform cell counts between 1 to 5 minutes after the addition of Trypan Blue to the cell suspension. Place 10 μ L of the cell suspension into at least one chamber of the hemocytometer. Count a minimum of 4 quadrants on the hemocytometer (see diagram below). Dead and dying cells are permeable to Trypan Blue and will look completely blue. For accurate cell counts, optimal number of cells per quadrant should be 25 to 75 cells. After counting the cells, calculate the average of the 4 quadrants. Take the cell count average and multiply by the dilution factor and by 10⁴ to get the number of cells per mL. Multiply the desired seeding density (5,000 viable cells per cm²) by the surface area of the vessel(s) to be inoculated. This will give you the total number of cells in the cell suspension. This will give you the volume of cells uspension with which to inoculate each vessel. Inoculate the cells into the culture vessels prepared with pre-warmed culture medium. Mix gently to evenly distribute the cells and place culture vessels into the incubator at 37°C, 5% CO₂.

Sample calculation:

Average viable cells per quadrant = 31

31 cells/quadrant x 10,000 quadrants/mL x 2 (dilution factor*) = 620,000 cells/mL *if using equal volumes of Trypan Blue and cell suspension.

Inoculating a T-75 flask at 5,000 cells/cm²:

5,000 cells/cm² x 75 cm² = 375,000 cells/flask

Calculate volume of cell suspension to inoculate each flask with:

375,000 cells/flask divided by 620,000 cells/mL = 0.604 mL/flask

Passaging Cells

Human Airway Epithelial Cells may be passaged once the culture is 70 to 80% confluent and actively proliferating. Human Airway Epithelial Cells are contact inhibited, and it is essential that the cells be passaged before reaching confluence, since post-confluent cells exhibit slower proliferation and morphological changes after passaging.

Lifeline® recommends using the TrypKit[™] Subculture Reagent Kit (LL-0013). Aspirate the medium from the culture vessel. Rinse the flasks with Lifeline's Phosphate Buffered Saline (PBS; CM-0001) by adding at least 1 mL of PBS per each 5 cm² and gently tilting the flask to cover the surface with PBS. Aspirate the PBS from the culture vessel, repeat the rinse if desired. Add at least 2 mL of Lifeline's 0.05% Trypsin/0.02% EDTA (CM-0017) per each 25 cm² to the vessel. Swirl gently to ensure all cells are coated with the Trypsin/EDTA. Observe the cells carefully under the microscope. When the cells round up they are ready to be released. This normally takes from 2 to 3 minutes depending on the confluence of the cells. Do not over trypsinize as this may damage the cells. Detach the cells by gently striking the culture vessel against your hand several times. Observe the cells under the microscope to be sure they have become detached. Once the cells are fully detached, add Lifeline's Trypsin Neutralizing Solution (TNS;



Call Lifeline Technical Service and Sales at 877.845.7787 Or visit lifelinecelltech.com for more information

Lifeline Cell Technology | 8415 Progress Drive, Suite T | Frederick, MD 21701 ©2024 Lifeline Cell Technology, an International Stem Cell Company | All Rights Reserved



 FC-0016
 LL-0023

 FC-0106
 LM-0007

 FC-0035
 LS-1047

 FC-0011
 FC-0054

 FC-0065

CM-0018) using a volume equal to the amount of Trypsin/EDTA that was originally used. Gently swirl to ensure all of the trypsin solution is neutralized. Using aseptic laboratory techniques, pipette the cells into a sterile centrifuge tube. Collect the remaining cells by rinsing the culture vessel with at least 1 mL of PBS per each 5 cm² and pipetting the cells into the sterile centrifuge tube. Check culture vessel under the microscope for cells still attached and repeat steps if necessary to retrieve all the cells from the vessel. All steps must be completed under aseptic conditions in a biological safety cabinet.

Centrifuge the cells at $150 \times g^*$ for 3 to 5 minutes. For best results, calculate speed for individual centrifuge type. Time may also be centrifuge dependent. Do not over centrifuge cells as this will cause cell damage. After centrifugation, the cells should form a clean loose pellet. Please consult Lifeline's technical service department if issues arise from trypsinization or centrifugation.

Aspirate the neutralized trypsin solution from the centrifuge tube and resuspend the cell pellet in prewarmed BronchiaLife[™] Medium by gently pipetting up and down with a 2- or 5-mL pipette. Count cells using a hemocytometer, and re-plate at 5,000 cells per cm² (or use a split ratio of 1:4 to 1:10) in vessel containing the respective pre-warmed culture medium.

*To calculate RCF (x g)

RCF = 0.00001118 × (rpm)² × r r = rotational radius in centimeters rpm = rotations or revolutions per minute

Recommended Feeding Guidelines

Guidelines for a T-25 Flask. Adjust volumes according to culture surface area.

Every other day, remove medium and feed with 5 mL of fresh supplemented medium.

Most cultures which are 50% confluent will be ready for passage the following day and should be fed with 7 to 8 mL of supplemented medium.

<u>Do not use more than 10 mL of medium per 25 cm²</u> of culture surface to ensure that the media depth is at a level where gas diffusion will be sufficient to support the cells' requirements for oxygen.

The depth of the medium affects gas diffusion gradients through the culture medium to the cells. The volumes of medium recommended in this table result in a range of depths between 2 mm and 5 mm, which is compatible with general recommendations, 10 mL being at the maximum depth allowable (5 mm).



Call Lifeline Technical Service and Sales at 877.845.7787 Or visit lifelinecelltech.com for more information Lifeline Cell Technology | 8415 Progress Drive, Suite T | Frederick, MD 21701 ©2024 Lifeline Cell Technology, an International Stem Cell Company | All Rights Reserved